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
Long-read capture with Twist target enrichment system

Long-read capture with Twist target enrichment system
Tina Han, Holly Corbitt, Leonardo Arbizu, Esteban Toro, Chad Locklear.
Twist Bioscience, South San Francisco, California, USA

1. Introduction

Targeted sequencing aims to facilitate characterization of gene regions at a scale-per-read that is more accessible than whole genome sequencing. While short-read NGS sequencing has been shown to accurately and comprehensively interrogate complex genomic architecture, short-reads often cannot uniquely focus on single genes using PCR enrichment-based methods. Here we describe a method to leverage Twist Bioscience target enrichment workflow for gene panels sequenced with long-reads. This custom platform the panelist will be the design phase of a panel of genes, including both exons and introns. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns.

2. Methods

We designed gene panels of various target sizes, ranging from 0.2-20 kb. Panels were optimized using a proprietary algorithm to create balanced capture of exons and introns at the same time, reducing variation of off-target sequences. Our long-read hybrid capture protocol used 200-1000 bp fragments. Panels that were created using commercial platforms, including Illumina or Agilent, were used as a baseline for comparison. The panels were then tested using PCR to determine if the panels could be used for long-read capture. The panels were then tested using PCR to determine if the panels could be used for long-read capture. The panels were then tested using PCR to determine if the panels could be used for long-read capture.

3. Long-Read Pharmacogenomics (PGx) Panel

Pharmacogenomics, including HLA genes and CYP2D6, are notoriously difficult to genotype via amplicon or short-read sequencing as the high level of nucleotide heterogeneity and complex structure of these genes often leads to a common artifact of whole-genome sequencing and resequencing. We set forth to design a 20-gene pharmacogenomics panel to address the issues associated with the HLA region and CYP2D6. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns.

4. Long-Read Dark Genes Panel

The repetitive nature and polymorphic complexity of 389 medically relevant genes poses a challenge for their accurate analysis in clinical settings. But with long-read sequencing, we have demonstrated that the panel provides all gene coverage for 302 difficult-to-read genes, including those genes in “NGS dead zones” that are difficult to sequence or map with short-reads. These genes are referred to as short coverage or missing information, read cliffs, and read gaps. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns.

5. Conclusions

We demonstrate that our long-read capture method effectively enables comprehensive coverage of gene targets using a panel of genes with highly complex genomic architecture, including those genes in “NGS dead zones” that are difficult to sequence or map with short-reads. This long-read hybrid capture protocol can be adapted to a wide range of gene panels to address the challenges of short-read long-read sequencing. Distorted sequencing reads due to common artifacts of whole-genome sequencing and resequencing can be avoided by using long-read capture to address the issues associated with the HLA region and CYP2D6. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns. The Twist Bioscience panel will be the design phase of a panel of genes, including both exons and introns.

Panel Size	Number of Genes	Size of Genes (kb)	Mean Target per Panel	Mean Target Coverage	Depth of Coverage	Full Enrichment	Duplicate Rate	Mean-Read Length
22 Mb	389	4	0.95%	75%	38x (50x on 302)	84	3%	5.2 kb

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https://www.twistbioscience.com/ja/resources/poster/long-read-capture-twist-target-enrichment-system

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https://www.twistbioscience.com/ja/resources/webinar/long-read-capture-twist-target-enrichment